

# Digital Balance 2kg 0.01g

- 1) House: Robust ABS
- 2) Interface: Rs232C (connect with Printer or PC)
- 3) LCD With Backlight
- 4) More functions: Counting, Full Range Tare, Over-Load Protection, Zero Tracking, Automatic Turn-Off, Under Weighing, Low Power Indicator
- 5) Value stable time: 2 s
- 6) Power: 230/115Vac  $\pm$  15%



## Attachments

- English Manual
- Adapter
- Windshield (No)
- M1 Class Calibration Weight (Optional)
- Under Weighing Hook (Optional)
- Rechargeable battery (Optional)
- RS232C Interface

## Specification

Model	WA1002F	WA2002F	WA3002F	WA5002F	WA6002F	WA10002F	WA20002F
Capacity	0-100g	0-200g	0-300g	0-500g	0-600g	0-1000g	0-2000g
Readability	0.01g						
Repeatability	$\leq 0.01g$	$\leq 0.01g$	$\leq 0.01g$	$\leq 0.01g$	$\leq 0.01g$	$\leq 0.02g$	$\leq 0.02g$
Linearity	$\leq 0.02g$	$\leq 0.02g$	$\leq 0.02g$	$\leq 0.02g$	$\leq 0.02g$	$\leq 0.03g$	$\leq 0.03g$
Pan dim.	188×160 mm						
Display	LCD With Backlight						
Line Voltage	AC (230/115Vac $\pm$ 15%)/DC (Dry Battery Or Rechargeable Battery)						
N. W.	2kg						
G. W.	2.5kg						
Outline	270×190×75mm						
Packing	450×295×550mm/carton (Including 8 pcs)						

Branch: 14 Fawzy Soliman St., From El Haram St., Giza, Egypt P.O.Box: 127 Giza, Egypt Code:12515  
Tel. Fax:+(202) 33943004 / Mobile: +2 (01223415321-0111221096-01011811884)  
[www.photoscience-co.com](http://www.photoscience-co.com) / [info@photoscience-co.com](mailto:info@photoscience-co.com)