


# DIGITAL COLONY COUNTER

**Use :** After cultivating of bacteria such as colon bacilli, it counts the number of bacteria. It is very useful lab equipment with a microscope for physiological research, hospital clinical diagnosis.



Standard accessories	
Magnifier 5" 1ea	Stainless Needle point pen 1ea
Power cord 1ea	Making pen probe 1ea
100mm Dish holder 1ea	

## Features

- By using of button and a marking pen, 2 counter mode are possible.
- This unit main power is 220V.
- Easy to read 4 digits LED digital display
- Buzzer sounds every counting.
- A marking pen is usable for user s convenience.
- Wider view of 125 mm wide lens.
- Bacteria can be counted by picking with stainless needle point pen.

## Specification

<b>Model NO</b>	<b>PH-DCC560</b>
<b>Display type</b>	4 Digital LED Display
<b>Petri dish holder</b>	Accepts 100mm dish (removable)150mm (fixed)
<b>Zero reset button</b>	Yes
<b>Illumination</b>	Fluorescent lamp
<b>Dimensions(WxDxH,mm)</b>	310x290x100
<b>Mode of Counting</b>	Manual push button countingNeedle-point probingMarking pen probing
<b>Electric Requirements</b>	220VAC. 50/60Hz

Standard accessories
Magnifier 5" 1ea
Stainless Needle point pen 1ea
Making pen probe 1ea
Power cord 1ea

**Branch: 14 Fawzy Soliman St., From El Haram St., Giza, Egypt P.O.Box: 127 Giza, Egypt Code:12515**  
**Tel. Fax:+(202) 33943004 / Mobile: +2 (01223415321-0111221096-01011811884)**  
**www.photoscience-co.com / info@photoscience-co.com**