

MICRO CENTRIFUGE

MAX. RPM 6,000 RCF
. 2,000 xg

Simple and Compact designed - personal centrifuge

Electronic breaker to stop the machine when lid opens. Convenient and time-saving lid open/close.

Remarkably durable and long last.

Prevents samples from exterior particles inflow, Fixed angle rotor for 8 X 1.5/2,0 ml tubes,

PR strip rotor for 2 strips X 8 X 0.2 ml tube&

Available in clean benches and hoods & cold room



SPECIFICATION

Max. RPM / RCF	6000 rpm
Max. RCF	2,000 xg
Max. Capacity	1,5/2,0m1x8
Noise level	≤55dB
RPM Accuracy	±100
Power supply	100~220V / 50~60Hz
Dimension(W x D x H)	162 x154x 116 mm
Weight (without rotor)	1kg

Branch: 14 Fawzy Soliman St., From El Haram St., Giza, Egypt **P.O.Box:** 127 Giza, Egypt Code:12515

Tel. Fax:+(202) 33943004 / **Mobile:** +2 (01223415321-0111221096-01011811884)

www.photoscience-co.com / info@photoscience-co.com