

SB-120DTN Ultrasonic Cleaner



Features:

- Digital display and setting ultrasonic cleaning time
- Display working countdown time
- Display actual temperature of the cleaning tank
- Temperature setting range: Ambient temperature-80°C
- Timing range: 1-99 minutes adjustable
- Equipped with a special stainless steel basket, reduce sound cover
- Shell and cover are made by stainless steel
- Equipment controlled by Single-chip controller
- Last working parameters can be stored in memory after power-off

Main Technical Parameters:

Model	Inner tank	volume	Frequency	Power	Heating power	Adjustable temperature	Adjustable time	Drainage	Rack
	L/W/H (mm)	L	KHz	W	W	°C	min		√
SB-100DTN	300*150*100	4.5	40	100	120	Ambient-80	1-99	/	√
SB-120DTP	240*140*100	3	40	120	120	Ambient-80	1-99		√
SB-120DTN	240*140*100	3	40	120	120	Ambient-80	1-99	/	√
SB-120DTN	240*140*150	5	40	120	120	Ambient-80	1-99	/	√
SB-3200DTN	300*155*150	6	40	180	150	Ambient-80	1-99	/	√
SB-5200DTN	300*240*150	10	40	200	320	Ambient-80	1-99	√	√
SB-5200DTN	300*240*150	10	40	250	320	Ambient-80	1-99	/	√
SB-5200DTN	300*240*150	10	40	300	320	Ambient-80	1-99	√	√
SB-4200DTN	300*240*200	14.4	40	400	320	Ambient-80	1-99	√	√
SB25-12DTN	500*300*150	22.5	40	500	450	Ambient-80	1-99	√	√
SB25-12DTN	500*300*150	22.5	40	600	450	Ambient-80	1-99	√	√

<https://www.photoscience.co.com>